

# HiEndoXL™ Endothelial Cell Expansion Medium, Reduced Serum

**Product Code: AL517**

## Product description:

HiEndoXL™ Endothelial Cell Expansion Medium is a reduced serum medium used for *in vitro* cultivation and expansion of Human Umbilical Vein Endothelial Cells (HUVEC) and Human Umbilical Artery Endothelial cells (HUAEC). It contains basal medium (Part A) and endothelial cell growth supplement (Part B). Part A consists of inorganic, organic salts, amino acids, vitamins and sodium bicarbonate. Part B consists of growth factors and nutrients necessary for growth of large vessel endothelial cells. This medium and supplement is devoid of antibiotics and antimycotics.

## Products Required But Not Supplied

|  |        |
|--|--------|
| <b>1. Media Supplements</b>                  |        |
| Antibiotic-Antimycotic Solution 100X [or]    | A002   |
| Gentamicin-Amphotericin B solution 1000X     | A031   |
| <b>2. Reagents for Sub-culture</b>           |        |
| Dulbecco's Phosphate Buffered Saline (DPBS)  | TL1006 |
| Trypsin/EDTA Solution 1X                     | TCL007 |
| EnVzyme™ Easy                                | TCL137 |
| Trypan Blue 0.5% solution                    | TCL005 |
| Trypsin Inhibitor from Soyabean              | TCL068 |
| <b>3. Reagent for Coating Culture vessel</b> |        |
| 0.5% Gelatin solution in DPBS                | TCL109 |

## Directions:

- Thaw endothelial cell growth supplement (Part B) overnight at 2-8°C.  
Note: Precipitates in Part B after thawing are normal. Precipitates will not affect the performance of the medium.

- Disinfect the external surface of the bottles of part A and Part B by spraying with isopropyl alcohol before placing in a biosafety hood.
- Transfer the entire content of Part B to basal medium (Part A) under aseptic condition.  
Note: If desired, 5ml of antibiotic-antimycotic solution (A002) can be added to 500ml of complete medium.
- Tightly cap the bottle and swirl gently to ensure proper mixing.  
**Note:** Do not mix vigorously. Doing so will cause formation of foam.
- Store the complete medium at 2 - 8°C until use.

## Quality control:

### Appearance

Part A: Orangish red coloured clear solution  
Part B: Pale yellow coloured clear solution

### pH

7.00-7.60

### Osmolality in mOsm/Kg H<sub>2</sub>O

290.00-330.00

### Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

### Cultural Response

The medium is tested for optimal cell growth and proliferation of mesenchymal stem cells.

## Storage and shelf life:

Store basal medium at 2-8°C away from bright light. Store endothelial growth supplement at -20 °C. Use before expiry date given on the product label. Shelf life of the complete medium is 6 weeks.

Note: Freezing of the basal medium and complete medium is not recommended. Avoid repeated freezing and thawing of the growth supplement.

**Table 1 : Gelatin Coating of Culture Vessel**





|  |  |   | Time Required (approx.) |
|--|--|---|-------------------------|
|  |  | Key Points to Remember  |                         |
|  |  | <b>For uniform coating , make sure that the incubator is properly levelled</b>                                  |                         |
| Aseptically add 0.5% gelatin solution (TCL109)                 |  | <b>Refer Table 2 for recommended volumes of gelatin solution</b>  | 1 min                   |
| Incubate overnight at 37°C incubator                           |  |   | Overnight               |
| Aspirate gelatin solution with the help of pipette.            |  |   |                         |
| <b>GELATIN COATED CULTURE VESSEL IS READY FOR USE</b>          |  |   |                         |
| If vessel is not used immediately, store at 2-8° upto one week |  | <b>Flask should be kept with caps tightly closed and plates should be sealed with a parafilm during storage</b> |                         |

**Table 2 : Recommended Volumes of Gelatin Solution for Different Culture Vessels**

| Culture Vessel | Volume Per Well |
|----------------|-----------------|
| 96-well plate  | 75 µl           |
| 48-well plate  | 150 µl          |
| 24-well plate  | 300 µl          |
| 12-well plate  | 500 µl          |
| 6-well plate   | 1 ml            |
| T-25 Flask     | 5 ml            |
| T-75 Flask     | 10 ml           |





**Table 3 : Protocol for Thawing**

- Cryopreserved cells are supplied in liquid nitrogen dry vapour shipper (-150°C to -130°C).
- Upon receipt, immediately transfer the vial to the vapor phase of liquid nitrogen tank.
- Store it in the tank until further use. Cells must be processed at least in a BSL II hood.

|   |   | Key Points to Remember  | Time Required (approx.) |
|---|---|---|-------------------------|
| <b>1. Preparation of Culture Vessel</b>   |   |   |                         |
| a. Add 5ml of complete medium to a T-25 flask   |    | <b>Preparation of complete medium</b><br>AL517 (Part A 500 ml) + (Part B 20 ml) + A002 (5 ml)                               | 60 secs                 |
| b. Place the flask at 37°C to equilibrate the medium  |    |   | 30 mins                 |
| <b>2. Thawing Procedure</b>   |   | <b>Make sure water bath is set at 37°C before starting the thawing procedure</b>  |                         |
| a. Remove cryovial from the liquid nitrogen tank/ shipper wearing appropriate protective gear   |  | <b>Thawing should be AS FAST AS POSSIBLE to minimize cell damage</b>  |                         |
| b. Immediately thaw the vial partially by holding in a water bath at 37°C   |  | <b>DO NOT hold the vial in water bath for more than 90-120 secs</b><br><b>AVOID getting water up to the cap of the vial</b> | 90-120 secs             |
| c. Disinfect the vial by swabbing thoroughly with 70% isopropyl alcohol   |  |   | 10 secs                 |
| d. Add the cell suspension drop by drop to the T-25 flask containing the pre-warmed complete medium. Keep swirling the flask while adding the cell suspension |  | <b>Dropwise addition is required to prevent the cells from stress induced by exothermic reaction</b>                        | 30-60 secs              |




**Table 3 : Protocol for Thawing**

- Cryopreserved cells are supplied in liquid nitrogen dry vapour shipper (-150°C to -130°C).
- Upon receipt, immediately transfer the vial to the vapor phase of liquid nitrogen tank.
- Store it in the tank until further use. Cells must be processed at least in a BSL II hood.

|   |   | Key Points to Remember   | Time Required (approx.)    |
|---|---|--|----------------------------|
| e. Cap the flask and shake gently to ensure proper mixing and uniform distribution of cells in the medium |    |  | 10 secs                    |
| <b>3. Incubation</b>  |   |  |                            |
| a. Incubate the cells at 37°C and 5% CO <sub>2</sub>  |    | Check for cell attachment in 2-3 hrs   | 2-3 hrs                    |
| b. If more than 70-80% cells are attached, replace the medium with fresh medium                           |   | <b>Medium change after 2-3 hours is mandatory to remove traces of DMSO</b><br><br>If cells have not attached, centrifuge the cell suspension at 1000 rpm for 7-8 mins and resuspend in fresh medium                | 60-120 secs<br><br>7-8 min |
| c. Incubate the cells at 37°C and 5% CO <sub>2</sub>  |  |  | 3-5 days                   |
| <b>YOUR CELLS ARE READY TO SUB-CULTURE</b>  |   |  |                            |
| <b>4. Maintenance</b>   |   |  |                            |
| a. Monitor the cells every day  |   | Use the recommended freezing medium for cryopreservation of cells<br><b>Upto 50% Confluency:</b><br><b>Change the medium on alternate day</b><br><b>After 50% Confluency:</b><br><b>Change the medium everyday</b> |                            |
| b. Change the medium  |   |  |                            |
| c. Sub-culture once cells reach 70 - 80% confluence   |   |  |                            |



**Table 4 : Subculture**

- HUVEC/HUAEC can be sub-cultured at a seeding density of 5000-10,000 cells/cm<sup>2</sup>.
- Sub-culturing ratios can vary from 1:2 - 1:5
- A confluent T-25 flask of HUVEC/HUAEC yields approx.  $1.0 \times 10^6$  cells

|   |   | Key Points to Remember  | Time Required (approx.) |
|---|---|---|-------------------------|
| a. Aspirate entire medium and discard<br>DO NOT disturb the monolayer   |    |   | 60 secs                 |
| b. Wash the cells with 2-3 ml DPBS to remove residual medium<br>c. Aspirate off the DPBS and discard  |    | Prior to use, make sure that <b>dissociation solution is equilibrated to room temperature</b>   | 60 secs                 |
| d. Add 0.5 ml pre-warmed Trypsin-EDTA solution or 1mL prewarmed EnVzyme™ Easy solution<br>e. <b>Incubate the flask at 37°C</b><br><br>- Trypsin-EDTA solution dissociates HUVEC in approx. 30 sec<br><br>- EnVzyme™ Easy solution dissociates HUVEC in approx. 7-10 min |   | Gently rock the flask to ensure complete coverage of the dissociation solution over the cells<br><br><b>Exposing the cells to Trypsin for longer time leads to loss of cell viability</b><br><br>EnVzyme™ Easy is gentle on cells and longer exposure does not harm the cells. It does not require neutralization | 30 sec<br><br>7-10 mins |
| f. Microscopically monitor the flask<br>g. When the cells start rounding up, gently tap the flask to ensure complete detachment of cells  |   |   | 15 secs                 |
| h. When using Trypsin-EDTA, neutralize its action by adding equal amount of Soyabean Trypsin Inhibitor Solution (TCL068).<br>i. When using EnVzyme™ Easy, add 1mL complete medium<br>j. Pipette gently to get a homogenous mixture of cells                             |  | <b>Vigorous pipetting will stress the cells</b>   | 60 secs                 |

**Table 4 : Subculture**

- HUVEC/HUAEC can be sub-cultured at a seeding density of 5000-10,000 cells/cm<sup>2</sup>.
- Sub-culturing ratios can vary from 1:2 - 1:5
- A confluent T-25 flask of HUVEC/HUAEC yields approx.  $1.0 \times 10^6$  cells

|  |   | Key Points to Remember  | Time Required (approx.) |
|--|---|---|-------------------------|
| k. Count cells using hemocytometer<br>l. Seed at recommended seeding density in a new flask containing fresh complete medium. Refer to Table 5 |  | <b>DO NOT refrigerate cells after splitting</b><br><b>Seed immediately</b>  | 10-15 mins              |
| m. Incubate in a humidified incubator at 37°C and 5% CO <sub>2</sub>   |  |   | 48 hrs                  |
| <b>Maintenance</b>   |   |   |                         |
| a. Monitor the cells every day<br>b. Change the medium<br>c. Sub-culture once cells reach 70 - 80% confluence                                  |   | <b>Upto 50% Confluency:</b><br><b>Change the medium on alternate day</b><br><b>After 50% Confluency:</b><br><b>Change the medium everyday</b> |                         |

**Table 5 : Seeding Density**

| Flask | Recommended Seeding Density  | No. of Cells Per Flask | Volume of Medium (ml) |
|-------|------------------------------|------------------------|-----------------------|
| T-25  | 5000 cells/cm <sup>2</sup>   | $0.125 \times 10^6$    | 5 - 7                 |
|       | 10,000 cells/cm <sup>2</sup> | $0.25 \times 10^6$     | 5 - 7                 |

These are recommended seeding densities from literature and our studies. Higher seeding densities do not cause any harm to the cells and reduce the required population doublings per passage. Lower seeding densities may cause cells to lose viability, detach during culture and in general take more population doublings to reach confluence.

## Related products:

| Product name  | Code   | Packing   |
|---|--|---|
| HiFi™ Human Umbilical Vein Endothelial Cells (HUVEC)                      | CL002-0.5<br>CL002-T25   | 0.5 million cells/vial<br>1 T25cm <sup>2</sup> flask  |
| HiFi™ Human Umbilical Vein Endothelial Cells (HUVEC) Without BBE and VEGF | CL021-0.5<br>CL021-T25   | 0.5 million cells/vial<br>1 T25cm <sup>2</sup> flask  |
| Accutase™   | TCL075-1X100ML<br>TCL075-5X100ML<br>TCL075-1X500ML                                       | 1 x 100ml<br>5 x 100ml<br>1 x 500ml   |
| Trypsin-EDTA Solution 1X  | TCL033-5X100ML<br>TCL033-2X500ML<br>TCL033-6X500ML                                       | 5 x 100ml<br>2 x 500ml<br>6 x 500ml   |
| Trypsin Inhibitor from soybean 1X; Liquid                                 | TCL068-1X100ML<br>TCL068-5X100ML   | 1 x 100ml<br>5 x 100ml  |
| Dulbecco's Phosphate Buffered Saline                                      | TL1006-5X100ML<br>TL1006-2X500ML<br>TL1006-6X500ML<br>TL1006-18X500ML<br>TL1006-1X1000ML | 5 x 100ml<br>2 x 500ml<br>6 x 500ml<br>18 x 500ml<br>1 x 1000ml   |
| Antibiotic Antimycotic solution 100X, Liquid                              | A002-5X20ML<br>A002-5X50ML<br>A002-5X100ML   | 5 x 20ml<br>5 x 50ml<br>5 x 100ml   |
| Gentamycin Solution   | A005-5X20ML<br>A005-5X50ML   | 5 x 20ml<br>5 x 50ml  |
| HiFi™ Human Umbilical Artery Endothelial cells (HUAEC)                    | CL017-0.5<br>CL017-2X0.5<br>CL017-T25<br>CL017-2XT25                                     | 0.5 million cells/vial<br>2 x 0.5 million cells/vial<br>1 T25cm <sup>2</sup> flask<br>2 T25cm <sup>2</sup> flasks |

Revision: 1 / 2017

### Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.