

# HiMesoXL™ Mesenchymal Stem Cell Expansion Medium, Reduced Serum

**Product Code: AL519**

## Product description:

HiMesoXL™ Mesenchymal Stem Cell Expansion Medium, is designed for *in vitro* cultivation and expansion of Wharton's Jelly Mesenchymal Stem Cells (HWJ-MSC), Human Adipose Derived Mesenchymal Stem Cells (HAD- MSC) and Human Dental Pulp Stem Cells (HDP-SC) while maintaining them in an undifferentiated state. It contains basal medium (Part A) and growth supplement (Part B). Part A consists of inorganic, organic salts, amino acids, vitamins and sodium bicarbonate and is devoid of protein, hormones, antibiotics and antimycotics. Part B is reduced serum growth supplement necessary for growth of mesenchymal stem cells.

- Disinfect the external surface of the bottles of part A and Part B by spraying with isopropyl alcohol before placing in a biosafety hood.
- Transfer the entire content of Part B to basal medium (Part A) under aseptic condition.  
*Note: If desired, 5ml of antibiotic-antimycotic solution (A002) can be added to 500ml of complete medium.*
- Tightly cap the bottle and swirl gently to ensure proper mixing.  
*Note: Do not mix vigorously. Doing so will cause formation of foam.*
- Store the complete medium at 2 - 8°C until use. The complete medium can be directly used for culturing cells and no further addition of serum is required.

## Quality control:

### Appearance

Part A: Orangish red coloured clear solution  
Part B: Pale yellow coloured clear solution

### pH

7.00-7.60

### Osmolality in mOsm/Kg H<sub>2</sub>O

280.00-320.00

### Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

### Cultural Response

The medium is tested for optimal cell growth and proliferation of mesenchymal stem cells.

## Storage and shelf life:

Store basal medium at 2-8°C away from bright light. Store mesenchymal cell growth supplement at -20°C.

Use before expiry date given on the product label. Shelf life of the complete medium is 6 weeks at 2-8°C.

**Note:** Do not freeze the basal medium. Avoid repeated freezing and thawing of the growth supplement.

## Products Required But Not Supplied

1.Media Supplements	Code
Mesenchymal Stem Cell Tested	RM10832, RM10845
Fetal Bovine Serum (FBS)	RM10846, RM10938
Antibiotic-Antimycotic Solution 100X[or] Gentamicin-Amphotericin B solution 1000X	A002 A031
2.Reagents for Sub-culture	Code
Dulbecco's Phosphate Buffered Saline (DPBS)	TL1006
Trypsin/EDTA Solution 1X	TCL007
Trypan Blue 0.5% solution	TCL005
Trypsin inhibitor from soyabean	TCL068
3. Stem Cell Freezing Medium	Code
CryoXL™ Stem Cell Freezing Medium	TCL107







## Directions:

- Thaw mesenchymal cell growth supplement (Part B) overnight at 2-8°C.

*Note: Few precipitates in Part B may be visible after thawing. Precipitates will not affect the performance of the medium.*





**Table 1 : Protocol for Thawing**

- Cryopreserved cells are supplied in liquid nitrogen dry vapour shipper (-150°C to -130°C).
- Upon receipt, immediately transfer the vial to the vapor phase of liquid nitrogen tank.
- Store it in the tank until further use. Cells must be processed at least in a BSL II hood.

		Key Points to Remember	Time Required (approx.)
<b>1. Preparation of Culture Vessel</b>			
a. Add 5ml of complete medium to a T-25 flask		<b>Preparation of complete medium</b> AL519 (Part A 500 ml) + (Part B 11.4 ml) + A002 (5 ml)	60 secs
b. Place the flask at 37°C to equilibrate the medium			30 mins
<b>2. Thawing Procedure</b>		Make sure <b>water bath</b> is set at <b>37°C</b> before starting the thawing procedure	
a. Remove cryovial from the liquid nitrogen tank/ shipper wearing appropriate protective gear		Thawing should be <b>AS FAST AS POSSIBLE</b> to minimize cell damage	
b. Immediately thaw the vial partially by holding in a water bath at 37°C		<b>DO NOT hold the vial</b> in water bath for more than 90-120 secs <b>AVOID</b> getting water upto the cap of the vial	90-120 secs
c. Disinfect the vial by swabbing thoroughly with 70% isopropyl alcohol			10 secs
d. Add the cell suspension drop by drop to the T-25 flask containing the pre-warmed complete medium. Keep swirling the flask while adding the cell suspension		<b>Dropwise addition</b> is required to <b>prevent the cells from stress</b> induced by exothermic reaction	30-60 secs





**Table 1 : Protocol for Thawing**

- Cryopreserved cells are supplied in liquid nitrogen dry vapour shipper (-150°C to -130°C).
- Upon receipt, immediately transfer the vial to the vapor phase of liquid nitrogen tank.
- Store it in the tank until further use. Cells must be processed at least in a BSL II hood.



		Key Points to Remember	Time Required (approx.)
e. Cap the flask and shake gently to ensure proper mixing and uniform distribution of cells in the medium			10 secs
<b>3. Incubation</b>			
a. Incubate the cells at 37°C and 5% CO <sub>2</sub>		Check for cell attachment in 2-3 hrs	2-3 hrs
b. If more than 70-80% cells are attached, replace the medium with fresh medium		<b>Medium change after 2-3 hours is mandatory to remove traces of DMSO</b>  If cells have not attached, centrifuge the cell suspension at 1000 rpm for 7-8 mins and resuspend in fresh medium	60-120 secs  7-8 min
c. Incubate the cells at 37°C and 5% CO <sub>2</sub>			3-5 days
<b>YOUR CELLS ARE READY TO SUB-CULTURE</b>			
<b>4. Maintenance</b>			
a. Monitor the cells every day b. Change the medium every alternate day c. Sub-culture, once cells reach 70 - 80% confluence		Use the recommended freezing medium for cryopreservation of cells <b>DO NOT</b> allow cells to reach 100% confluency before sub culture or cryopreservation  In case of reduced serum or serum free media , use trypsin inhibitor solution (TCL068) for neutralisation of Trypsin during subculture  Usage of just medium for neutralisation will result in inefficient neutralisation and will stress the cells resulting in reduced viability and cell death	

**Table 2 : Subculture**

- HWJ-MSC /HAD-MSC/HDP-SC can be sub-cultured at a seeding density of 5000-10,000 cells/cm<sup>2</sup>.
- Sub-culturing ratios can vary from 1:2 - 1:5
- A confluent T-25 flask of HWJ-MSC /HAD-MSC/HDP-SC yields 1.0 x 10<sup>6</sup> cells

		Key Points to Remember	Time Required (approx.)
a. Aspirate entire medium and discard DO NOT disturb the monolayer			60 secs
b. Wash the cells with 2-3 ml DPBS to remove residual medium c. Aspirate off the DPBS and discard		Prior to use, make sure that <b>Trypsin-EDTA solution is equilibrated to room temperature</b>	60 secs
d. Add 0.5 ml pre-warmed Trypsin-EDTA solution		Gently rock the flask to ensure complete coverage of the Trypsin-EDTA solution over the cells	
e. Incubate the flask in the incubator at 37°C for 30-60 secs		<b>Exposing the cells to Trypsin-EDTA for longer time leads to loss of cell viability</b>	30-60 secs
f. Microscopically monitor the flask. g. When the cells start rounding up, gently tap the flask to ensure complete detachment of cells			15 secs
h. To neutralize action of trypsin add 3 ml of complete medium, if AL512 is used i. Pipette gently to get a homogenous mixture of cells j. If reduced serum medium AL519 is used, add 0.5 ml Soyabean Trypsin Inhibitor Solution (TCL068). Centrifuge the cell suspension at 1000 rpm for 10 mins. Discard supernatant and resuspend pellet in fresh 3 ml of complete medium by pipetting.		<b>Vigorous pipetting will stress the cells</b>	60 secs

## Table 2 : Subculture

<ul style="list-style-type: none"> <li>• HWJ-MSC /HAD-MSC/HDP-SC can be sub-cultured at a seeding density of 5000-10,000 cells/cm<sup>2</sup>.</li> <li>• Sub-culturing ratios can vary from 1:2 - 1:5</li> <li>• A confluent T-25 flask of HWJ-MSC /HAD-MSC/HDP-SC yields 1.0 x 10<sup>6</sup> cells</li> </ul>			
		Key Points to Remember	Time Required (approx.)
<p>k. Count cells using hemocytometer</p> <p>l. Seed at recommended seeding density in a new flask containing fresh complete medium</p> <p>Refer to Table 3</p>		<p><b>DO NOT refrigerate</b> cells after splitting</p> <p><b>Seed immediately</b></p>	10-15 mins
<p>m. Incubate in a humidified incubator at 37°C and 5% CO<sub>2</sub></p>			48 hrs
<b>Maintenance</b>			
<p>a. Monitor the cells every day</p> <p>b. Change the medium every alternate day</p> <p>c. Sub-culture once cells reach 70 - 80% confluence</p>			

## Table 3 : Seeding Density

Flask	Recommended Seeding Density	No. of Cells Per Flask	Volume of Medium (ml)
T-25	5000 cells/cm <sup>2</sup>	0.125 x 10 <sup>6</sup>	5 - 7
	10,000 cells/cm <sup>2</sup>	0.25 x 10 <sup>6</sup>	5 - 7

These are recommended seeding densities from literature and our studies. Higher seeding densities do not cause any harm to the cells and reduce the required population doublings per passage. Lower seeding densities may cause cells to lose viability, detach during culture and in general take more population doublings to reach confluence.

**Related products:**

<b>Product Name</b>	<b>Code</b>	<b>Packing</b>
HiFi™ Wharton's Jelly Derived Mesenchymal Stem cells (HWJ-MSC)	CL001-0.5	0.5 million cells/vial
	CL001-T25	1 T25cm <sup>2</sup> flask
	CL001-T75	1 T75cm <sup>2</sup> flask
HiFi™ Adipose Derived Mesenchymal Stem cells (HAD-MSC)	CL007-0.5	0.5 million cells/vial
	CL007-T25	1 T25cm <sup>2</sup> flask
	CL007-T75	1 T75cm <sup>2</sup> flask
HiFi™ Human Dental Pulp Stem cells (H-DPSC)	CL008-0.5	0.5 million cells/vial
	CL008-T25	1 T25cm <sup>2</sup> flask
	CL008-T75	1 T75cm <sup>2</sup> flask
HiMesoXL™ Mesenchymal Stem Cell Expansion Medium, Reduced serum	AL519-1X100ML	1x100ml
	AL519-1X500ML	1x500ml
HiAdipoXL™ Adipocyte Differentiation Medium	AL521-1X100ML	1x100ml
HiOsteoXL™ Osteocyte Differentiation Medium	AL522-1X100ML	1x100ml
HiChondroXL™ Chondrocyte Differentiation Medium	AL523-1X100ML	1x100ml
EZXpand™ Mesenchymal Stem Cell Culture Kit (Adipose-derived)	CCK024-0.5	0.5 million cells/vial
	CCK024-T25	1 T25cm <sup>2</sup> flask
EZXpand™ Mesenchymal Stem Cell Culture Kit (Wharton's Jelly derived)	CCK025-0.5	0.5 million cells/vial
	CCK025-T25	1 T25cm <sup>2</sup> flask
CryoXL™ Stem Cell Freezing Medium	TCL107-1X50ML	1x50ml
Accutase™	TCL075-1X100ML	1 x 100ml
	TCL075-5X100ML	5 x 100ml
	TCL075-1X500ML	1 x 500ml
Trypsin-EDTA Solution 1X	TCL033-5X100ML	5 x 100ml
	TCL033-2X500ML	2 x 500ml
	TCL033-6X500ML	6 x 500ml
Trypsin Inhibitor from soybean 1X; Liquid	TCL068-1X100ML	1 x 100ml
	TCL068-5X100ML	5 x 100ml
Dulbecco's Phosphate Buffered Saline	TL1006-5X100ML	5 x 100ml
	TL1006-2X500ML	2 x 500ml
	TL1006-6X500ML	6 x 500ml
	TL1006-18X500ML	18 x 500ml
	TL1006-1X1000ML	1 x 1000ml
Antibiotic Antimycotic solution 100X, Liquid	A002-5X20ML	5 x 20ml
	A002-5X50ML	5 x 50ml
	A002-5X100ML	5 x 100ml
Gentamicin Solution w/ 50mg/ml Gentamicin in sterile tissue culture grade water	A005-5X20ML	5 x 20ml
	A005-5X50ML	5 x 50ml

Revision: 2 / 2017

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.