

PCT0611 Maltodextrin

Product Number Packing

PCT0611 : 500G

Product Information

Product Code : PCT0611
Product Name : Maltodextrin
Shelf Life : 4 years

Technical Specification

Appearance : White to almost white hygroscopic granules or powder
Solubility : 100 mg soluble in 1 mL of water
Cultural response : Cultures conditions - Incubation period (5wks), Relative humidity (60±2%), Temperature (25±2°C), Photoperiod Day: Night in hours (16:8)
Shoot culture : No structural deformity observed, actively growing shoots, no toxicity to shoots
Callus culture : No necrotic tissues, actively growing callus, no toxicity to callus
Identification : Complies (To 5 mL of hot alkaline cupric tartrate TS add a few drops of a 1 in 10 solution of the sample. A red precipitate of cuprous oxide form)
Reducing sugar : ≤ 20.00%

Risk and Safety Information

WGK : 3
Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : -
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods