

Antibiotic Antimycotic Solution 100X liquid,

Endotoxin tested
With 10,000 U Penicillin, 10mg Streptomycin and 25 µg Amphotericin B per ml in 0.9% normal saline

Product Code: A002A

Product Description :

Penicillin and streptomycin are broad spectrum bacteriostatic and bacteriocidal agents effective against Gram positive and Gram negative bacteria. Amphotericin B is an antifungal agent effective against fungi and yeasts. Penicillin acts by inhibiting cell wall synthesis. Streptomycin blocks the initiation complex of protein synthesis by binding to the 30S subunit of bacterial 70S ribosome, thus interfering with the protein synthesis. Amphotericin B interferes with the cell membrane permeability of sensitive fungi and yeasts by forming channels and causing leakage of low molecular weight substances.

A002A is a sterile filtered solution formulated to contain 10,000 units penicillin, 10mg streptomycin and 25µg amphotericin B per ml.

Directions :

Recommended use concentration is 10ml/L. Sensitive cell lines may react differently to this product. Hence it is recommended to determine optimum usage dose empirically for individual cell line.

Quality Control:

Appearance

Yellow suspension.

pH

4.00 to 7.00

Osmolality in mOsm/Kg H₂O

325.00 -365.00

Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

Penicillin Concentration

10,000 - 12,000 U/ml.

Streptomycin Concentration

10 -12 mg/ml.

Amphotericin Concentration

25 - 30 ug/ml.

Cultural Response

1. No toxicity to cells.
2. Antibiotic sensitivity test

Endotoxin Content

NMT 0.5EU/ml

Storage and Shelf Life:

Store at -20°C. Repeated freezing and thawing should be avoided. Once thawed, remaining portion can be aseptically dispensed into sterile container and frozen for future use.

Solution is stable for 3 days at 37°C and for approximately 2 to 3 weeks at 2 - 8°C.

Shelf life of the product is 24 months.

Note: Precipitate in this product is normal. Slight precipitate may be visible if the bottles are not adequately warmed. However, the precipitate gets fully dissolved again on warming the solution at room temperature.

Use before expiry date given on the product label.

Disclaimer:

Revision: 4/ 2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.