www.himedialabs.com

Safety data sheet(SDS)

According to Regulation (EC) No.1907/2006

Revision : 01

Date of Revision : 18.10.2017

1 Identification of the substances/ mixture and of the company/ undertaking

1.1	Product Identifiers						
	Product Number	A018					
	Product Name	Antibiotic Solution 100X liquid					
		w/ 10,000U Penicillin, 5mg Streptomycin per ml in 0.9% normal saline					
1.2	Relevant identified uses of	the substance or mixture and uses advise	d against				
1.2.1	Relevant identified uses	Laboratory chemicals, Manufacture of sub	ostances				
1.3	Details of the supplier of the safety data sheet						
	Produced by	by HiMedia Laboratories Private Limited					
	Address	23, Vadhani Industrial Estate, LBS Marg, Ghatkopar (W), Mumbai - 400 086 India					
	Tel. No.	+91-22-2500 0970, +91-22-2500 1607	Fax No. : +91-22-25002468				
	Mail Id	info@himedialabs.com	Website : www.himedialabs.com				
1.4	Emergency Tel. No.						
	Emergency Tel. No.	Please contact the regional HiMedia repr	esentation in your country				

2 Hazards Identification

HIMEDIA

2.1 Classification of the substance or mixture *CLP Classification-Regulation (EC) No. 1272/2008[EU-GHS/CLP]*

Not a hazardous substance or mixture according to Regulation (EC) No.1272/2008

2.2 Label elements Labeling according to Regulation (EC) No.1272/2008

The product does not need to be labelled in accordance with EC directives or respective national laws.

2.3 Other Hazards

None

3 Composition/Information On Ingredients

No components need to be disclosed according to the applicable regulations.

4 First Aid Measures

4.1 Description of first aid measures

General advice

Show this safety data sheet to the doctor in attendance.

If inhaled

If breathed in, move person into fresh air. Consult a physician.

In case of skin contact

Wash off with soap and plenty of water. If skin irritation occurs,get medical advice/attention.

Page **1** of **6**

In case of eye contact Rinse out with plenty of water with the eyelid held wide open. If eye irritation persists, get medical advice/attention. If swallowed Rinse mouth with water. Consult a physician if feeling unwell. 4.2 Most important symptoms and effects, both acute and delayed The most important known symptoms and effects are described in the labeling (see section 2.2) and/or in section 11. 4.3 Indication of immediate medical attention and special treatment needed No data available 5 **Fire Fighting Measures** 5.1 **Extinguishing media** Suitable extinguishing media Use extinguishing measures that are appropriate to local circumstances and the surrounding environment. Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide. Unsuitable extinguishing media 5.2 Special hazards arising from the substance or mixture No data available 5.3 **Precautions for fire-fighters** No data available. 5.4 **Further information** Wear self-contained breathing apparatus for firefighting if necessary. 6 **Accidental Release Measures** 6.1 Personal precautions, protective equipment and emergency procedures Use personnel protective equipment. Wear disposable gloves, dust mask and eye protection. Avoid dust formation. For personal protection see section 8. 6.2 **Environmental precautions** Prevent further leakage or spillage if safe to do so. No special environmental precautions required. 6.3 Methods and materials for containment and cleaning up Keep in suitable, closed containers for disposal. 6.4 **Reference to other sections** For disposal see Section 13. Handling and Storage 7 7.1 Precautions for safe handling Avoid contact with skin and eyes. Avoid formation of dust and aerosols. Provide appropriate exhaust ventilation at places where dust is formed. For precautions see section 2.2. 7.2 Conditions for safe storage, including any incompatibilities Store in cool place. Keep container tightly closed in a dry and well-ventilated place. Containers which are opened must be carefully resealed and kept upright to prevent leakage. Recommended Storage Temperature : store between -5 to -20° C Page 2 of 6

7.3 Specific end uses

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated.

8 Exposure Controls/Personal Protection

8.1 Control parameters

Components with workplace control parameters

8.2 Exposure controls

Appropriate engineering controls

Handle in accordance to general industrial hygiene and safety practice. Wash hands before breaks and immediately after handling the products.

Personal protective equipment

Hygiene measure

Avoid contact with skin, eyes and clothing. Immediately change contaminated clothing.

Eye/face protection

Safety glasses with side-shields conforming to EN 166 Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166 (EU).

Skin protection

The selected protective gloves have to satisfy the specifications of EU Directive 89/686/EEC and the standard EN 374 derived from it. Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Body protection

Impervious clothing The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

Respiratory protection

For nuisance exposures use type P95 (US) or type P1 (EU EN 143) particle respirator. For higher level protection use type OV/AG/P99 (US) or type ABEK-P2 (EU EN 143) respirator cartridges. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Environment exposure controls

Do not let product enter drains.

9 Physical and chemical properties

9.1 Information on basic physical and chemical properties

Appearance Odour Odour Threshold pH Melting/freezing point Initial boiling point and boiling range Flash point Upper/lower flammability or explosive limits Evaporation rate Flammability (Solid, gas) Vapour pressure Relative density

Colorless, clear solution No data available Partition coefficient: n-octanol/water Autoignition Temperature **Decomposition Temperature** Viscosity **Explosive properties Oxidizing properties** Vapour density Thermal decomposition

No data available No data available

9.2 Other safety information

No data available

10	Stability and Reactivity
10.1	Reactivity
	No data available
10.2	Chemical stability
	Stable under recommended storage conditions.
10.3	Possibility of hazardous reactions
	No data available
10.4	Conditions to avoid
	No data available
10.5	Incompatible materials
	No data available
10.6	Hazardous decomposition products
	None under normal use conditions. Other decomposition products. No data available. In event of fire
	- refer section 5

11

No data available

Toxicological Information 11.1 Information on toxicological effects Acute toxicity No data available Skin corrosion/irritation No data available Serious eye damage/eye irritation No data available Respiratory or skin sensitisation No data available Germ cell mutagenicity No data available Carcinogenicity IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC. Reproductive toxicity No data available

Specific target organ toxicity- single exposure

Page 4 of 6

Specific target organ toxicity - repeated exposure No data available Aspiration hazard No data available

12 **Ecological Information**

12.1 Toxicity

No data available

- 12.2 Persistence and degradability No data available
- 12.3 **Bioaccumulative potential** No data available
- 12.4 Mobility in soil
- No data available
- 12.5 PBT and vPvB assessment PBT and vPVB assessment not available Chemical safety assessment is not required
- 12.6 Other adverse effects No data available

13 **Disposal Considerations**

13.1 Waste treatments methods

Product

Offer surplus and non-recyclable solutions to a licenced disposal company. Dispose off waste in accordance with all applicable Federal, state and local laws.

13.2 Contaminated packaging

Dispose in accordance with all applicable federal, state, and local environmental regulations.

14 **Transport Information**

14.1 UN-No

ADNR : ADR : IATA_C : IATA_P : IMDG : RID :

14.2 UN proper shipping name

- : Not dangerous goods ADNR
- ADR : Not dangerous goods
- : Not dangerous goods
- IATA_C IATA_P : Not dangerous goods
- IMDG : Not dangerous goods
- RID : Not dangerous goods

14.3 Transport hazard class(es)

ADNR : - ADR : - IATA_C : - IATA_P : - IMDG : - RID : -

14.4 Packaging group

Page 5 of 6

	ADNR : ADR :	IATA_C :	IATA_P :	IMDG :	RID :	
14.5 14.6	Environmental hazards Special precautions for use No data available					
15	Regulatory Information This safety datasheet compl	ies with the require	ments of Regulat	ion(EC) No. 1907/	2006.	
15.1	mixture					
	mixture No data available					

16 Other information

Further Information

Copyright 2016 HiMedia Laboratories Pvt. Ltd.

The information given in this safety data sheet is believed to be correct yet does not claim to be all inclusive. This document is intended only as a guide for appropriate precautionary handling of the material by properly trained individuals, information here being commensurate with the present state of our knowledge regarding the manner and conditions of use, handling, storage or disposal. The information provided herein shall not be considered as guarantee of the properties of the product. HiMedia Laboratories, shall not be held liable for any damage resulting from improper handling or contact with the above product. Unless explicitly stated on the product or in any of the documentation accompanying the product, it is intended for research or testing purpose only and is not to be used for any other purpose.