



HiGlutaXL[™] Dulbecco's Modified Eagle Medium, High Glucose

With 4.5gms Glucose per litre, L-Alanyl-L-Glutamine, Sodium pyruvate and Sodium bicarbonate

Product Code: AL007G

Product Description:

HiGlutaXLTM medium contains the stabilised dipeptide form of L-glutamine, L-alanyl-L-glutamine. HiGlutaXLTM medium offers several advantages over the conventional glutamine containing media. Dipeptide form prevents the intramolecular cyclization reaction, thus preventing toxic build -up of ammonia. L-alanyl-Lglutamine incorporates L-alanine that protects the alpha amino acid group. Aminopeptidases within the cell break the dipeptide, gradually releasing both L-glutamine and L-alanine for use by the cell. The gradual release of Lglutamine obviates the need to supplement L-glutamine frequently and allows liquid media to be stored at 4°C for longer periods.

AL007G is HiGlutaXLTM Dulbecco's Modified Eagle Medium with 4.5gms glucose per liter, L-alanyl-L-glutamine, sodium pyruvate and sodium bicarbonate. Users are advised to review the literature for recommendations regarding medium supplementation and physiological growth requirements specific for different cell lines.

Composition:

Ingredients	mg/L
INORGANIC SALTS	
Calcium chloride dihydrate	265.000
Ferric nitrate nonahydrate	0.100
Magnesium sulphate anhydrous	97.720
Potassium chloride	400.000
Sodium bicarbonate	3700.000
Sodium chloride	6400.000
Sodium dihydrogen phosphate anhydrous	109.000
AMINO ACIDS	
Glycine	30.000
L-Alanyl-L-Glutamine	868.000
L-Arginine hydrochloride	84.000
L-Cystine dihydrochloride	62.570
L-Histidine hydrochloride monohydrate	42.000
L-Isoleucine	105.000
L-Leucine	105.000

L-Lysine hydrochloride	146.000
L-Methionine	30.000
L-Phenylalanine	66.000
L-Serine	42.000
L-Threonine	95.000
L-Tryptophan	16.000
L-Tyrosine disodium salt dihydrate	103.790
L-Valine	94.000
VITAMINS	
Choline chloride	4.000
D-Ca-Pantothenate	4.000
Folic acid	4.000
Nicotinamide	4.000
Pyridoxal hydrochloride	4.000
Riboflavin	0.400
Thiamine hydrochloride	4.000
i-Inositol	7.200
OTHERS	
D-Glucose	4500.000
Phenol red sodium salt	15.900
Sodium pyruvate	110.000

Quality Control:

Appearance

Red colored, clear solution.

pН

7.00 -7.60

Osmolality in mOsm/Kg H₂O

320.00 - 360.00

Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

Cultural Response

The growth promotion capacity of the medium is assessed qualitatively by analyzing the cells for the morphology and quantitatively by estimating the cell counts and comparing it with a control medium.

Endotoxin Content

NMT 1EU/ml

Storage and Shelf Life:

Store at 2-8 °C away from bright light. Shelf life is 18 months. Use before expiry date given on the product label.

Disclaimer: Revision: 02/2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

