



10.000

20.000

60.000

70.000

15.000

25.000

# **Technical Datasheet**

# Medium 199

# With Earle's salts, L-Glutamine and Sodium bicarbonate

**Product Code: AL014A** 

# **Product Description:**

Medium 199 was the first nutritionally defined medium developed by Morgan, Morton, and Parker in 1950. This complex medium was formulated specifically for nutritional studies on primary chick embryo fibroblasts in the absence of any additives. It was observed that explanted tissue could survive in Medium 199 without serum but long term cultivation of cells required of the medium supplementation with serum. Medium 199 is formulated with either Hank's salts or Earle's salts. The medium when supplemented with serum can be used for growth of a wide variety of cells. Medium 199 is presently used for the maintenance of non-transformed cells, vaccine and virus production and explants epithelial primary of cells. AL014A is Medium 199 with Earle's salts, L-glutamine and sodium bicarbonate. Users are advised to review the literature for recommendations regarding medium supplementation and physiological growth requirements specific for different cell lines.

## **Composition:**

Ingredients	mg/L
INORGANIC SALTS	2 < 7 000
Calcium chloride dihydrate	265.000
Magnesium sulphate anhydrous	97.720
Potassium chloride	400.000
Sodium acetate anhydrous	50.000
Sodium bicarbonate	2200.000
Sodium chloride	6800.000
Sodium phosphate monobasic	122.000
AMINO ACIDS	0.050
Ascorbic acid	0.050
Glycine	50.000
L-Alanine	25.000
L-Arginine hydrochloride	70.000
L-Aspartic acid	30.000
L-Cysteine hydrochloride monohydrate	0.100
L-Cystine dihydrochloride	26.000
L-Glutamic acid	67.000
L-Glutamine	100.000
L-Histidine hydrochloride monohydrate	22.000

<b>,</b>	
L-Proline	40.000
L-Serine	25.000
L-Threonine	30.000
L-Tryptophan	10.000
L-Tyrosine disodium salt dihydrate	57.660
L-Valine	25.000
VITAMINS	0.400
Calciferol	0.100
Choline chloride	0.500
D-Biotin	0.010
D-Ca-Pantothenate	0.010
DL-Tocopherol phosphate disodium salt	0.010
Folic acid	0.010
Menadione	0.010
Nicotinamide	0.025
Nicotinic acid	0.025
Pyridoxal hydrochloride	0.025
Pyridoxine hydrochloride	0.025
Retinol Acetate	0.140
Riboflavin	0.010
Thiamine hydrochloride	0.010
i-Inositol	0.050
p-Amino benzoic acid (PABA)	0.050
OTHERS	40.000
Adenine sulphate	10.000
Adenosine monophosphate	0.200
Adenosine triphosphate	1.000
Cholesterol	0.200
Deoxyribose	0.500
Ferric nitrate nonahydrate	0.720
Glucose	1000.000
Glutathione reduced	0.050
Guanine hydrochloride	0.300
Hypoxanthine	0.354
Phenol red	15.000

L-Hydroxyproline

L-Lysine hydrochloride

L-Isoleucine

L-Methionine

L-Phenylalanine

L-Leucine

Polysorbate	4.900
80 Ribose	0.500
Thymine	0.300
Uracil	0.300
Xanthine	0.344

## **Quality Control:**

## Appearance

Orangish red colored, clear solution.

#### рH

7.00 - 7.60

## Osmolality in mOsm/Kg H<sub>2</sub>O

280.00 -320.00

### **Sterility**

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

## **Cultural Response**

The growth promotion capacity of the medium is assessed qualitatively by analyzing the cells for the morphology and quantitatively by estimating the cell counts.

## **Endotoxin Content**

NMT 1EU/ml

## **Storage and Shelf Life:**

Store at 2-8°C away from bright light. Shelf life is 12 months. Use before expiry date given on the product label.

Disclaimer: Revision: 04/2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic , research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

