



HiGlutaXL™ Medium 199

With L-Alanyl-L-Glutamine, Earle's Salts and Sodium bicarbonate
Without HEPES buffer

Product Code: AL014G

Product description:

HiGlutaXL™ medium contains the stabilised dipeptide form of L-glutamine, L-alanyl-L-glutamine. HiGlutaXL™ medium offers several advantages over the conventional glutamine containing media. Dipeptide form prevents the intramolecular cyclization reaction, thus preventing toxic build-up of ammonia. L-alanyl-L-glutamine incorporates L-alanine that protects the alpha amino acid group. Aminopeptidases within the cell break the dipeptide, gradually releasing both L-glutamine and L-alanine for use by the cell. The gradual release of L-glutamine obviates the need to supplement L-glutamine frequently and allows liquid media to be stored at 4°C for longer periods.

AL014G is HiGlutaXL™ Medium 199 with L-alanyl-L-Glutamine, Earle's salts and sodium bicarbonate. Users are advised to review the literature for recommendations regarding medium supplementation and physiological growth requirements specific for different cell lines.

Composition:

Ingredients	mg/L
INORGANIC SALTS	
Calcium chloride dihydrate	265.000
Ferric nitrate nonahydrate	0.720
Magnesium sulphate anhydrous	97.720
Potassium chloride	400.000
Sodium acetate anhydrous	50.000
Sodium bicarbonate	2200.000
Sodium chloride	6800.00
Sodium phosphate monobasic	122.00
AMINO ACIDS	
Glycine	50.000
L-Alanine	25.000
L-Alanyl-L-Glutamine	148.782
L-Arginine hydrochloride	70.000
L-Aspartic acid	30.000
L-Cysteine hydrochloride monohydrate	0.100
L-Cystine dihydrochloride	26.00

L-Glutamic acid	67.000
L-Histidine hydrochloride monohydrate	22.000
Hydroxyproline	10.000
L-Isoleucine	20.000
L-Leucine	60.000
L-Lysine hydrochloride	70.000
L-Methionine	15.000
L-Phenylalanine	25.000
L-Proline	40.000
L-Serine	25.000
L-Threonine	30.000
L-Tryptophan	10.000
L-Tyrosine disodium salt dihydrate	57.660
L-Valine	25.000

VITAMINS:	
Ascorbic acid	0.050
Calciferol	0.100
Choline chloride	0.500
D-Biotin	0.010
D-Ca-Pantothenate	0.010
DL-Tocopherol phosphate disodium salt	0.010
Folic acid	0.010
Menadione	0.010
Nicotinamide	0.025
Nicotinic acid	0.025
Pyridoxal hydrochloride	0.025
Pyridoxine hydrochloride	0.025
Retinol acetate	0.140
Riboflavin	0.010
Thiamine hydrochloride	0.010
i-Inositol	0.050
p-Amino benzoic acid (PABA)	0.050
OTHERS	
Adenine sulphate	10.000
Adenosine monophosphate	0.200
Adenosine triphosphate	1.000
Cholesterol	0.200

Deoxyribose	0.5000
Glucose	1000.000
Glutathione reduced	0.050
Guanine hydrochloride	0.300
Hypoxanthine sodium salt	0.354
Phenol red sodium salt	15.000
Polysorbate 80	4.900
Ribose	0.500
Thymine	0.300
Uracil	0.300
Xanthine	0.344

Quality control:

Appearance

Orangish red coloured clear solution

pH

7.00 – 7.60

Osmolality in mOsm/Kg H₂O

280.00 – 320.00

Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification

Cultural Response

The growth promotion capacity of the medium is assessed qualitatively by analyzing the cells for the morphology and quantitatively by estimating the cell counts.

Endotoxin Content

NMT 1EU/ml

Storage and shelf life:

Store at 2-8°C away from bright light.

Shelf life is 18 months.

Use before expiry date given on the product label.

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Revision: 03/2022