



Technical Datasheet

Minimum Essential Medium Eagle (MEM)

With Hanks' salts, L-Glutamine and Sodium bicarbonate

Product Code: AL056A

Product Description:

Minimum Essential Medium (MEM) is a modification of Basal Medium Eagle (BME). It was developed by Harry Eagle to meet the specific nutritional requirements of certain subtypes of HeLa cells and normal mammalian fibroblasts. MEM includes higher concentration of amino acids so as to closely approximate the protein composition of cultured mammalian cells. MEM can be used either with Earle's salts or Hanks' salts and can also be additionally supplemented with Non-essential Amino Acids (NEAA). This medium can be further modified by eliminating calcium to facilitate growth of cells in suspension cultures.

AL056A is Minimum Essential Medium with Hanks' balanced salts, L-glutamine and sodium bicarbonate. Hanks' salt mixture is designed to equilibrate with air, hence does not require CO₂ air mixture. Cells can therefore be grown in AL056A in less CO₂ or CO₂ free environment. Users are advised to review the literature for recommendations regarding medium supplementation and physiological growth requirements specific for different cell lines.

Composition:

| Ingredients | mg/L |
|---------------------------------------|----------|
| INORGANIC SALTS | |
| Calcium chloride dihydrate | 185.000 |
| Disodium hydrogen phosphate anhydrous | 47.800 |
| Magnesium sulphate anhydrous | 97.720 |
| Potassium chloride | 400.000 |
| Potassium dihydrogen phosphate | 60.000 |
| Sodium bicarbonate | 350.000 |
| Sodium chloride | 8000.000 |
| AMINO ACIDS | |
| L-Arginine hydrochloride | 126.000 |
| L-Cystine dihydrochloride | 31.300 |
| L-Glutamine | 292.000 |
| L-Histidine hydrochloride monohydrate | 42.000 |
| L-Isoleucine | 52.000 |
| | |

| L-Leucine | 52.000 |
|------------------------------------|----------|
| L-Lysine hydrochloride | 72.500 |
| L-Methionine | 15.000 |
| L-Phenylalanine | 32.000 |
| L-Threonine | 48.000 |
| L-Tryptophan | 10.000 |
| L-Tyrosine disodium salt dihydrate | 51.900 |
| L-Valine | 46.000 |
| VITAMINS | |
| Choline chloride | 1.000 |
| D-Ca-Pantothenate | 1.000 |
| Folic acid | 1.000 |
| Nicotinamide | 1.000 |
| Pyridoxal hydrochloride | 1.000 |
| Riboflavin | 0.100 |
| Thiamine hydrochloride | 1.000 |
| i-Inositol | 2.000 |
| OTHERS | |
| D-Glucose | 1000.000 |
| Phenol red sodium salt | 11.000 |
| | |

Quality Control:

Appearance

Red colored, clear solution

pН

7.00 - 7.60

Osmolality in mOsm/Kg H₂O

310.00 - 350.00

Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

Cultural Response

The growth promotion capacity of the medium is assessed qualitatively by analyzing the cells for the morphology and quantitatively by estimating the cell counts.

Endotoxin Content

NMT 1EU/ml

Storage and Shelf Life:

Store at 2-8°C away from bright light. Shelf life is 12 months. Use before expiry date given on the product label.

Disclaimer: Revision: 04/2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.