



Technical Datasheet

Waymouth Medium MB 752/1

With L-Glutamine and Sodium bicarbonate

Product Code: AL091A

Product Description:

Waymouth Medium MB 752/1 was developed in 1959 for the cultivation of mouse L cells, NCTC clone 929. Waymouth's medium is a comparatively simple medium with lesser ingredients than its other synthetic counterparts. The medium comprises of amino acids, vitamins, the purine base hypoxanthine, salts, and glucose. This medium can be used for culturing many fastidious cell lines and also whole organ culture. AL091A is Waymouth Medium MB 752/1 with L-glutamine and sodium bicarbonate. Users are advised to review the literature for recommendations regarding medium supplementation and physiological growth requirements specific for different cell lines.

Composition:

Ingredients	mg/L
INORGANIC SALTS	
Calcium chloride dihydrate	120.000
Disodium hydrogen phosphate anhydrous	300.000
Magnesium chloride anhydrous	112.560
Magnesium sulphate anhydrous	97.720
Potassium chloride	150.000
Potassium dihydrogen phosphate	80.000
Sodium bicarbonate	2240.000
Sodium chloride	6000.000
AMINO ACIDS	
Glycine	
L-Arginine hydrochloride	50.000
L-Aspartic acid	
L-Cysteine hydrochloride	75.000
L-Cystine dihydrochloride monohydrate	60.000 90.000
L-Glutamic acid	19.550
L-Glutamine	
L-Histidine hydrochloride monohydrate	150.000
L-Isoleucine	350.000
L-Leucine	150.000
	25.000
L-Lysine hydrochloride	50.000
L-Methionine	240.000
L-Phenylalanine	50.000
L-Proline	50.000
	50.000

L-Threonine	75.000
L-Tryptophan	40.000
L-Tyrosine disodium salt dihydrate	57.660
L-Valine	65.000
VITAMINS	
Ascorbic acid	17.500
Biotin	0.020
Choline chloride	250.000
D-Ca-Pantothenate	1.000
Folic acid	0.400
Niacinamide	1.000
Pyridoxine hydrochloride	1.000
Riboflavin	1.000
Thiamine hydrochloride	10.000
Vitamin B12	0.200
i-Inositol	1.000
OTHERS	
D-Glucose	5000.000
Glutathione reduced	15.000
Hypoxanthine sodium salt	29.000
Phenol red sodium salt	11.000

Quality Control:

Appearance

Red colored, clear solution

pН

7.00 - 7.60

Osmolality in mOsm/Kg H2O

300.00 -340.00

Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

Cultural Response

The growth promotion capacity of the medium is assessed qualitatively by analyzing the cells for the morphology and quantitatively by estimating the cell counts.

Endotoxin Content

NMT 1EU/ml



Storage and Shelf Life:

Store at 2-8°C away from bright light. Shelf life is 12 months. Use before expiry date given on the product label.

Disclaimer: Revision: 04/2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

