



Technical Datasheet

HiGlutaXL[™] Dulbecco's Modified Eagle Medium/Nutrient Mixture F-12 Ham (DMEM/ F12, 1:1 mixture)

With L-Alanyl-L-Glutamine, Sodium bicarbonate and Trace elements Without HEPES buffer

Product Code: AL127G

Product Description:

HiGlutaXLTM medium contains the stabilised dipeptide form of L-glutamine, L-alanyl-L-glutamine. HiGlutaXLTM medium offers several advantages over the conventional glutamine containing media. Dipeptide form prevents the intramolecular cyclization reaction, thus preventing toxic build -up of ammonia. L-alanyl-Lglutamine incorporates L-alanine that protects the alpha amino acid group. Aminopeptidases within the cell break the dipeptide, gradually releasing both L-glutamine and L-alanine for use by the cell. The gradual release of Lglutamine obviates the need to supplement L-glutamine frequently and allows liquid media to be stored at 4°C for longer periods.

AL127G is HiGlutaXLTM DMEM/Nutrient Mixture F-12 Ham with L-alanyl-L-glutamine, sodium bicarbonate and trace elements. It does not contain HEPES buffer. Users are advised to review the literature for recommendations regarding medium supplementation and physiological growth requirements specific for different cell lines.

Composition:

Ingredients	mg/L
INORGANIC SALTS	-
Ammonium metavanadate	0.00058
Ammonium molybdate tetrahydrate	0.00618
Calcium chloride dihydrate	154.500
Copper sulphate pentahydrate	0.0013
Disodium hydrogen phosphate	71.020
Ferric nitrate nonahydrate	0.050
Ferrous sulphate heptahydrate	0.417
Magnesium chloride hexahydrate	61.200
Magnesium sulphate anhydrous	48.840
Manganese sulphate	0.000151
Nickel chloride	0.00012
Potassium chloride	311.800
Sodium bicarbonate	1200.000

Sodium chloride	6996.000
Sodium dihydrogen phosphate	54.300
Sodium metasillicate nonahydrate	0.0142
Sodium selenite	0.00519
Stannous chloride dihydrate	0.00011
Zinc sulphate heptahydrate	0.432
AMINO ACIDS	
Glycine	18.750
L-Alanine	4.450
L-Alanyl-L-Glutamine	543.000
L-Arginine hydrochloride	147.500
L-Asparagine monohydrate	7.500
L-Aspartic acid	6.650
L-Cysteine hydrochloride monohydrate	17.560
L-Cystine dihydrochloride	31.290
L-Glutamic acid	7.350
L-Histidine hydrochloride monohydrate	31.480
L-Isoleucine	54.470
L-Leucine	59.050
L-Lysine hydrochloride	91.250
L-Methionine	17.240
L-Phenylalanine	35.480
L-Proline	17.250
L-Serine	26.250
L-Threonine	53.450
L-Tryptophan	9.020
L-Tyrosine disodium salt dihydrate	48.100
L-Valine	52.850
VITAMINS	0.000
Choline chloride	8.980
D-Biotin	0.0035
D-Ca-Pantothenate	2.240
Folic acid	2.660
Niacinamide	2.020
Pyridoxal hydrochloride	2.000
Pyridoxine hydrochloride	0.031

Riboflavin

0.219

Thiamine hydrochloride Vitamin B12 myo-Inositol	2.170 0.680 12.600
OTHERS	3151.000
D-Glucose	3131.000
DL-Thioctic acid	0.105
Hypoxanthine sodium salt	2.400
Linoleic acid	0.042
Phenol red sodium salt	8.630
Putrescine hydrochloride	0.081
Sodium pyruvate	110.000
Thymidine	0.365

Quality Control:

Appearance

Red colored, clear solution.

pН

7.00 - 7.60

Osmolality in mOsm/Kg H₂O

290.00 -330.00

Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

Cultural Response

The growth promotion capacity of the medium is assessed qualitatively by analyzing the cells for the morphology and quantitatively by estimating the cell counts.

Endotoxin Content

NMT 1EU/ml

Storage and Shelf Life:

Store at 2-8°C away from bright light. Shelf life is 18 months.

Use before expiry date given on the product label.

Disclaimer: Revision: 03/2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

