



Technical Datasheet

Nutrient Mixture F-12 Ham, Coon's Modification

With Sodium bicarbonate Without L-Glutamine

Product Code: AL146

Product Description:

Ham's Nutrient Mixtures were originally developed for single cell plating of near diploid Chinese hamster ovary (CHO) cells and mouse L-cells. Both F-10 and F-12 are formulated for use with or without serum, depending on the type of cells being cultured.

Ham's Nutrient Mixture F12 was originally designed for serial propagation and cloning of two CHO cell lines namely, CHD-3 and CHL-1 and mouse L cells. It is the medium of choice for the growth of cells of rodent origin and for cloning of myeloma and hybridoma cells. This medium is also the medium of choice for clonal toxicity assay using CHO cells.

Ham's Nutrient Mixture F12, Coon's modification was developed to determine optimal conditions for obtaining viable hybrids by treatment of parental cells with sendai virus. The medium was modified by doubling the concentration of amino acids and pyruvate and adding ascorbic acid.

AL146 is Coon's modification of Nutrient Mixture F-12 Ham without L-glutamine.

Users are advised to review the literature for recommendations regarding medium supplementation and physiological growth requirements specific for different cell lines.

Composition:

Ingredients	mg/L
INORGANIC SALTS	
Calcium chloride dihydrate	165.000
Copper sulphate pentahydrate	0.0025
Disodium hydrogen phosphate, anhydrous	132.040
Ferric sulphate heptahydrate	0.834
Magnesium chloride hexahydrate	106.000
Magnesium sulphate anhydrous	50.800
Potassium chloride	305.000
Potassium dihydrogen phosphate	68.000
Sodium bicarbonate	1176.000
Sodium chloride	7530.000
Zinc sulphate heptahydrate	0.144

AMINO ACIDS

Sodium pyruvate

Thymidine

Glycine	16.000
L-Alanine	18.000
L-Arginine hydrochloride	422.000
L-Asparagine monohydrate	30.000
L-Aspartic acid	26.000
L-Cysteine hydrochloride monohydrate	70.260
L-Glutamic acid	30.000
L-Histidine hydrochloride monohydrate	42.000
L-Isoleucine	7.800
L-Leucine	26.200
L-Lysine hydrochloride	73.000
L-Methionine	9.000
L-Phenylalanine	10.000
L-Proline	70.000
L-Serine	21.000
L-Threonine	23.800
L-Tryptophan	4.000
L-Tyrosine disodium salt dihydrate	15.860
L-Valine	23.400
VITAMINS	
Biotin	0.0073
Choline chloride	13.960
D-Ca-Pantothenate	0.480
Folic acid	1.320
Nicotinamide	0.04
Pyridoxine hydrochloride	0.06
Riboflavin	0.04
Thiamine hydrochloride	0.340
Vitamin B12	1.360
i-Inositol	18.020
OTHERS	
D-Glucose	2000.000
Hypoxanthine sodium salt	4.040
Linoleic acid	0.09
Lipoic acid	0.206
Phenol red sodium salt	1.240
Putrescine dihydrochloride	0.300

220.100

0.700

Directions:

1. Add 10ml of 200mM L-glutamine (TCL012) for 1 litre of medium.

Material required but not provided:

L-Glutamine solution 200mM (TCL012)

Quality Control:

Appearance

Light orange colored, clear solution

nН

7.00 - 7.60

Osmolality in mOsm/Kg H₂O

320.00 - 360.00

Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

Cultural Response

The growth promotion capacity of the medium is assessed qualitatively by analyzing the cells for the morphology and quantitatively by estimating the cell counts.

Endotoxin Content

NMT 1EU/ml

Storage and Shelf Life:

Store at 2-8°C away from bright light. Shelf life is 18 months. Use before expiry date given on the product label.

Disclaimer: Revision: 04/2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic , research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

