



# **RPMI-1640**

With L-Glutamine, 2gms per liter Glucose, 0.165 moles per liter MOPS buffer and Sodium bicarbonate

# **Product Code: AL180A**

# **Product Description:**

Roswell Park Memorial Institute (RPMI) media are a series of media developed by Moore et al for the culture of human normal and neoplastic cells in vitro. RPMI 1640 is the most commonly used medium in the series. A modification of McCoy's 5A medium, the medium was specifically designed to support the growth of human lymphoblastoid cells in suspension culture. Presently the medium is extensively used for a wide range of anchorage dependant cell lines. The medium needs to be supplemented with 5-20% fetal bovine serum. The medium is also known to support growth of cells in the absence of serum.

AL180A is RPMI 1640 with L-glutamine, 2gms per litre glucose, 0.165 moles per litre MOPS buffer and sodium bicarbonate. MOPS, a zwitterionic buffer does not antagonize antifungal agents at final concentration of 0.165mol/L for pH 7.0. Therefore, this medium is used as a diluent for antifungal agents that are water-soluble as well as waterinsoluble. For water-insoluble antifungal agents, that cannot be prepared as stock solutions in water, such as amphotericin B, anidulafungin, itraconazole, ketoconazole, posaconazole and voriconazole, a dilution series of the agent should be prepared first at 100 times final strength in an appropriate solvent. Each of these non-aqueous solutions should then be diluted tenfold **RPMI-1640** broth". in Users are advised to review the literature for recommendations regarding medium supplementation and physiological growth requirements specific for different cell lines. Adapted from Clinical and Laboratory Standards Institute document M27-A3 - Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Approved Standard - Third edition; Vol.28 No.14

### **Composition:**

Ingredients	mg/L
INORGANIC SALTS	
Calcium nitrate tetrahydrate	100.000

Magnesium sulphate anhydrous	48.840
Potassium chloride	400.000
Sodium bicarbonate	2000.000
Sodium chloride	6000.000
Sodium phosphate dibasic anhydrous	800.000
AMINO ACIDS	
Glycine	10.000
L-Arginine hydrochloride	241.000
L-Asparagine monohydrate	50.000
L-Aspartic acid	20.000
L-Cystine dihydrochloride	65.200
L-Glutamic acid	20.000
L-Glutamine	300.000
L-Histidine hydrochloride monohydrate	20.960
L-Hydroxyproline	20.000
L-Isoleucine	50.000
L-Leucine	50.000
L-Lysine hydrochloride	40.000
L-Methionine	15.000
L-Phenylalanine	15.000
L-Proline	20.000
L-Serine	30.000
L-Threonine	20.000
L-Tryptophan	5.000
L-Tyrosine disodium salt dihydrate	28.830
L-Valine	20.000
VITAMINS	2 000
Choline chloride	$3.000 \\ 0.200$
D-Biotin	0.250
D-Ca-Pantothenate	1.000
Folic acid	1.000
Niacinamide	1.000
Pyridoxine hydrochloride	0.200
Riboflavin	1.000
Thiamine hydrochloride	0.005
Vitamin B12	35.000
i-Inositol	1.000
p-Amino benzoic acid (PABA) OTHERS	1.000
D-Glucose	2000.000
Glutathione reduced	1.000
Staumone reduced	1.000

MOPS buffer, free acid	34500.000
Phenol red sodium salt	5.300

# **Quality Control:**

#### Appearance

Orangish red colored, clear solution

#### pН

7.00 -7.60

# Osmolality in mOsm/Kg H<sub>2</sub>O

470.00 -510.00

#### Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

#### **Cultural Response**

The growth promotion capacity of the medium is assessed qualitatively by analyzing the cells for the morphology and quantitatively by estimating the cell counts.

**Endotoxin Content** NMT 1EU/ml

## **Storage and Shelf Life:**

Store at 2-8°C away from bright light. Shelf life is 12 months. Use before expiry date given on the product label.

#### Disclaimer :

Revision : 04/2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia<sup>™</sup> publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia<sup>™</sup> Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic , research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

