

Leibovitz's L-15 Medium 2X

With L-Glutamine
1X Liquid Cell Culture Medium

Product Code: AL247A

Product Description :

Leibovitz's Medium was specifically designed to grow cells in a CO₂ free atmosphere. The standard sodium bicarbonate/ CO₂ buffering system is replaced by combination of free basic amino acids, phosphate buffers and higher levels of galactose and sodium pyruvate. As a result, the medium does not require supplementation with sodium bicarbonate and can be used under conditions of free gaseous exchange with the atmosphere. The medium can be used to grow human tumor cells and embryonic cells and also established cell lines like HeLa and Hep-2. The medium is frequently used in diagnostic virology where tissue cell lines or strains need to be grown in closed systems. Leibovitz's medium obviates the need of frequent medium change.

AL247A is 2X Leibovitz's L-15 Medium with L-glutamine. Users are advised to review the literature for recommendations regarding medium supplementation and physiological growth requirements specific for different cell lines.

Composition :

Ingredients	mg/L
INORGANIC SALTS	
Calcium chloride dihydrate	370.000
Magnesium chloride hexahydrate	400.000
Magnesium sulphate anhydrous	195.440
Potassium chloride	800.000
Potassium phosphate monobasic	120.000
Sodium chloride	16000.000
Sodium phosphate dibasic anhydrous	380.240
AMINO ACIDS	
DL-Alpha alanine	900.000
Glycine	400.000
L-Arginine (free base)	1000.000
L-Asparagine	500.000
L-Cysteine (free base)	240.000
L-Glutamine	600.000
L-Histidine (free base)	500.000

L-Isoleucine	500.000
L-Leucine	250.000
L-Lysine hydrochloride	188.000
L-Methionine	150.000
L-Phenylalanine	250.000
L-Serine	400.000
L-Threonine	600.000
L-Tryptophan	40.000
L-Tyrosine disodium salt	552.320
L-Valine	200.000
VITAMINS	
Choline chloride	2.000
D-Ca-Pantothenate	2.000
Folic acid	2.000
Nicotinamide	2.000
Pyridoxine hydrochloride	2.000
Riboflavin-5-phosphate sodium salt	0.200
Thiamine monophosphate	2.000
i-Inositol	4.000
OTHERS	
D-Galactose	1800.000
Phenol red sodium salt	22.000
Sodium pyruvate	1100.000

Quality Control:

Appearance

Orangish red coloured clear solution.

pH (at 2X concentration)

7.30 -7.90

Osmolality in mOsm/Kg H₂O (at 2X concentration)

620.00 -680.00

Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

Cultural Response

The growth promotion capacity of the medium is assessed qualitatively by analyzing the cells for the morphology and

quantitatively by estimating the cell counts and comparing it with a control medium.

Endotoxin Content

NMT 1EU/ml

Storage and Shelf Life:

Store at 2-8°C away from bright light.

Shelf life is 12 months.

Use before expiry date given on the product label.

Revision : 0 / 2016

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.