



Nutrient Mixture F-12 Ham, Kaighn's Modification

With 1.5 gms per litre Sodium bicarbonate Without L-Glutamine

Product Code: AL271S

Product Description:

Ham's Nutrient Mixtures were originally developed for single cell plating of near diploid Chinese hamster ovary (CHO) cells and mouse L-cells. Both F-10 and F-12 are formulated for use with or without serum, depending on of the type cells being cultured. Ham's Nutrient Mixture F12 was originally designed for serial propagation and cloning of two CHO cell lines namely, CHD-3 and CHL-1 and mouse L cells. It is the medium of choice for the growth of cells of rodent origin and for cloning of myeloma and hybridoma cells. This medium is also the medium of choice for clonal toxicity CHO assay using cells. Kaighn's modification of Ham's F-12 is a complex formulation of F-12 with increased amounts of amino acids and puruvate. Salts used in this formulation are as given by Konisberg.

AL271S is Nutrient Mixture F-12 Ham with sodium bicarbonate. It does not contain L-glutamine. Users are advised to review the literature for recommendations regarding medium supplementation and physiological growth requirements specific for different cell lines.

Composition:

Ingredients	mg/L
INORGANIC SALTS	135.240
Calcium chloride dihydrate	
Copper sulphate pentahydrate	0.0025
Disodium hydrogen phosphate	115.020
Ferrous sulphate heptahydrate	0.834
Magnesium chloride anhydrous	49.700
Magnesium sulphate anhydrous	93.700
Potassium chloride	283.290
Potassium dihydrogen phosphate	58.500
Sodium chloride	7597.200
Sodium bicarbonate	1500.00
Zinc sulphate hepthydrate	0.1437

AMINO ACIDS	
Glycine	15.010
L-Alanine	17.800
L-Arginine hydrochloride	421.400
L-Asparagine monohydrate	30.020
L-Aspartic acid	26.620
L-Cystine hydrochloride monohydrate	70.240
L-Glutamic acid	29.420
L-Histidine hydrochloride monohydrate	41.920
L-Isoleucine	7.872
L-Leucine	26.240
L-Lysine hydrochloride	73.040
L-Methionine	8.960
L-Phenylalanine	9.920
L-Proline	69.060
L-Serine	21.020
L-Threonine	23.820
L-Tryptophan	4.080
L-Tyrosine disodium salt dihydrate	13.500
L-Valine	23.420
VITAMINS	
Biotin	0.073
Choline chloride	13.960
D-Ca-Pantothenate	0.477
Folic acid	1.320
Niacinamide	0.037
Pyridoxine hydrochloride	0.061
Riboflavin	0.0376
Thiamine hydrochloride	0.337
Vitamin B12	1.355
i-Inositol	18.020
OTHERS	
D-Glucose	1260.000
Hypoxanthine sodium salt	4.083
Lipoic acid	0.2063
Phenol red Sodium Salt	3.318
Putrescine dihydrochloride	0.322
Sodium pyruvate	220.000
Thymidine	0.726

Directions:

1. Add 10ml of 200mM L-glutamine (TCL012) for 1 litre of medium.

Material required but not provided:

L-Glutamine solution 200mM (TCL012)

Quality Control:

Appearance

Light orange colored, clear solution

pН

7.00 -7.60

Osmolality in mOsm/Kg H₂O 320.00 -360.00

Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

Cultural Response

The growth promotion capacity of the medium is assessed qualitatively by analyzing the cells for the morphology and quantitatively by estimating the cell counts.

Endotoxin Content NMT 1EU/ml

Storage and Shelf Life:

Store at 2-8°C away from bright light. Shelf life is 18 months. Use before expiry date given on the product label.

Disclaimer :

Revision: 04/2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic , research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

