



# HiEndoXL<sup>™</sup> Microvascular Endothelial Cell Expansion Medium, Reduced Serum

Without Bovine Brain Extract

## Product Code: AL527

## **Product description:**

HiEndoXLTM Microvascular Endothelial Cell Expansion Medium is a reduced serum medium used for in vitro cultivation and expansion of Human Microvascular Endothelial Cells. It contains Basal Medium (Part A) and Microvascular Endothelial Growth Supplement (Part B). Part A consists of inorganic, organic salts, amino acids, vitamins and sodium bicarbonate.Part B consists of growth factors and nutrients necessary for growth of large vessel endothelial cells. This medium and supplement is devoid of Bovine Brain Extract (BBE), antibiotics and antimycotics.

## **Products Required But Not Supplied**

1Media Supplements	Code	
Mesenchymal Stem Cell Tested	RM10832, RM10845	
Fetal Bovine Serum (FBS)	RM10846, RM10938	
Antibiotic-Antimycotic Solution 100X[or] Gentamicin-	A002	
Amphotericin B solution 1000X	A031	
2.Reagents for Sub-culture	Code	
Dulbecco's Phosphate	TL1006	
Buffered Saline (DPBS)		
Trypsin/EDTA Solution 1X	TCL007	
Trypan Blue 0.5% solution	TCL005	
Trypsin inhibitor from soyabean	vabean TCL068	
3. Stem Cell Freezing Medium	Code	
CryoXL™ Stem Cell Freezing	TCL107	
Medium		

## **Directions:**

1. Thaw melanocyte growth supplement (Part B) overnight at 2-8°C.

*Note: Precipitates in Part B after thawing are normal. Precipitates will not affect the performance of the medium.* 

- 2. Disinfect the external surface of the bottles of part A and Part B by spraying with isopropyl alcohol before placing in a biosafety hood.
- 3. Transfer the entire content of Part B to basal medium (Part A) under aseptic condition.

*Note: If desired, 5ml of antibiotic-antimycotic solution* (A002) can be added to 500ml of complete medium.

4. Tightly cap the bottle and swirl gently to ensure proper mixing.

**Note**: Do not mix vigorously. Doing so will cause formation of foam.

5. Store the complete medium at 2 - 8°C until use.

## **Quality control:**

#### Appearance

Part A: Orangish red coloured clear solution Part B: Pale yellow coloured clear solution

### pН

7.00 - 7.60

Osmolality in mOsm/Kg H2O 280.00 -320.00

#### Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

#### **Cultural Response**

The medium is tested for optimal cell growth and proliferation of microvascular endothelial cells.

### Storage and shelf life:

Store basal medium at 2-8°C away from bright light. Store endothelial growth supplement at -20 °C. Use before expiry date given on the product label. Shelf life of the complete medium is 6 weeks at 2-8°C. **Note:** Freezing of the basal medium and complete medium is not recommended. Avoid repeated freezing and thawing of the growth supplement.

Product Name	Code	Packing
HiEndoXL <sup>TM</sup> Endothelial Growth	TCL112-1NO	For 500ml of basal
Supplement, Reduced Serum w/o BBE	102112 1110	medium
HiEndoXL <sup>TM</sup> Large Vessel Endothelial Growth Supplement, Reduced Serum w/o BBE	TCL113-1NO	For 500ml of basal medium
HiFi <sup>TM</sup> Umbilical Vein Endothelial Cells(HUVEC	CL002-0.5 CL002-T25 CL002-T75	0.5 million cells/vial 1 T25cm <sup>2</sup> flask 1 T75cm <sup>2</sup> flask
EZXpand <sup>™</sup> Umbilical Vein Endothelial Cell Culture Kit	CCK026-0.5 CCK026-T25	0.5 million cells/vial 1 T25cm <sup>2</sup> flask
CryoXL <sup>TM</sup> Stem Cell Freezing Medium	TCL107-1X50ML	1x50ml
Accutase <sup>TM</sup>	TCL075-1X100ML TCL075-5X100ML	1 x 100ml 5 x 100ml
	TCL075-1X500ML	1 x 500ml
Trypsin-EDTA Solution 1X	TCL033-5X100ML TCL033-2X500ML	5 x 100ml 2 x 500ml
Trypsin-EDTA Solution TA	TCL033-6X500ML	6 x 500ml
Trypsin Inhibitor from soybean 1X; Liquid	TCL068-1X100ML TCL068-5X100ML	1 x 100ml 5 x 100ml
Dulbecco's Phosphate Buffered Saline	TL1006-5X100ML TL1006-2X500ML TL1006-6X500ML TL1006-18X500ML TL1006-1X1000ML	5 x 100ml 2 x 500ml 6 x 500ml 18 x 500ml 1 x 1000ml
Antibiotic Antimycotic solution 100X, Liquid	A002-5X20ML A002-5X50ML A002-5X100ML	5 x 20ml 5 x 50ml 5 x 100ml
Gentamicin Solution w/ 50mg/ml Gentamicin in sterile tissue culture grade water	A005-5X20ML A005-5X50ML	5 x 20ml 5 x 50ml
Gentamicin – Amphotericin B Solution 1000X	A031-1X20ML A031-5X20ML	1 x 20ml 5 x 20ml

Revision: 1 / 2015

#### **Disclaimer:**

**Related products:** 

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia<sup>™</sup> Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia<sup>™</sup> Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.