



Basal Medium Eagle (BME)

(Modified for Autoclaving) With Earle's salts Without L-Glutamine and Sodium bicarbonate

Product Code: AT004A

Product Description :

Basal Medium developed by Harry Eagle is a combination of essential nutrients in appropriate concentrations for monolayer cultivation of a wide variety of normal and transformed cells. The medium was initially developed as a result of studies to determine the nutritional requirements of HeLa cells and mouse fibroblast L cells in culture. Although there are many versions of Basal Medium described by Eagle, the name Basal Medium Eagle applies to only the formulation developed for HeLa cells. Basal Medium Eagle when properly supplemented supports growth of variety of diploid or primary mammalian cell cultures. Modifications to the original BME have resulted in other media, including MEM and DMEM.

AT004A is Basal Medium Eagle with Earle's balanced salts and is modified for autoclaving. Autoclavable media offer a convenient alternative to membrane sterilized liquid medium. It is modified to include heat stable components to ensure that product efficacy is maintained after autoclaving. L-glutamine is heat labile, hence has been omitted from the formulation.

Users are advised to review the literature for recommendations regarding medium supplementation and physiological growth requirements specific for different cell lines.

Composition:

Ingredients	mg/L
INORGANIC SALTS	
Calcium chloride dihydrate	265.000
Magnesium sulphate anhydrous	97.720
Potassium chloride	400.000
Sodium chloride	6800.000
Sodium dihydrogen phosphate anhydrous	122.000
AMINO ACIDS	
L-Arginine hydrochloride	21.100
L-Cystine dihydrochloride	15.650
L-Histidine hydrochloride monohydrate	10.500
L-Isoleucine	26.200

L-Leucine L-Lysine hydrochloride L- Methionine L-Phenylalanine L-Threonine L-Tryptophan L-Tyrosine disodium salt dihyrate L-Valine VITAMINS Choline bitartarate D-Biotin D-Ca-Pantothenate Folic acid	$\begin{array}{c} 26.200\\ 36.480\\ 7.500\\ 16.500\\ 23.800\\ 4.000\\ 25.950\\ 23.400\\ 1.800\\ 1.000\\ 1.000\\ 1.000\\ 1.000\\ 1.000\\ 1.000\\ \end{array}$
Nicotinamide Pyridoxal hydrochloride Riboflavin Thiamine hydrochloride i-Inositol OTHERS D-Glucose Phenol red disodium salt Sodium succinate Succinic acid	$\begin{array}{c} 1.000\\ 0.100\\ 1.000\\ 1.800\\ \end{array}$ $\begin{array}{c} 1000.000\\ 6.360\\ 100.000\\ 75.000\\ \end{array}$

Directions :

1. Suspend 9.1gms in 900 ml tissue culture grade water with constant, gentle stirring until the powder is completely dissolved. Do not heat the water.

2. Adjust the pH to 4.0 before autoclaving.

3. Make up the volume to 960ml (This volume is derived after subtracting the volumes of 7.5% sodium bicarbonate solution and 200mM L-glutamine solution from the final volume).

4. Autoclave medium at 121°C at 15psi for 15 minutes.

5. Remove the medium promptly from the autoclave to avoid extended heating or evaporation.

6. Allow to cool at room temperature.

7. Add 29.3ml of 7.5% sodium bicarbonate solution (TCL013) and 10 ml of 200mM L- glutamine solution (TCL012) to the final volume of the medium being prepared. 8. If necessary, adjust the pH using sterile 1N NaOH (TCL002) or 1N HCl (TCL003).

9. Store liquid medium at 2-8°C and in dark till use.

Material required but not provided :

Tissue culture grade water (TCL010) Sodium bicarbonate (TC230) Sodium bicarbonate solution, 7.5% (TCL013) L-Glutamine (TC243) 200mM L-Glutamine solution (TCL012) 1N Hydrochloric acid (TCL003) 1N Sodium hydroxide (TCL002) Foetal bovine serum (RM1112/RM10432)

Quality Control:

Appearance

White to light pink, homogenous powder

Solubility

Clear solution at 9.1 gms/L.

pH without Sodium Bicarbonate 4.30 -4.90

pH with Sodium Bicarbonate 7.20 -7.80

Osmolality without Sodium Bicarbonate(mOsm/Kg H₂O) 230.00 -270.00

Osmolality with Sodium Bicarbonate(mOsm/Kg H₂O) 280.00 -320.00

Cultural Response

The growth promotion capacity of the medium is assessed qualitatively by analyzing the cells for the morphology and quantitatively by estimating the cell counts.

Endotoxin Content NMT 1EU/ml

Storage and Shelf Life:

1. All the powdered media and prepared liquid culture media should be stored at 2-8°C. Use before the expiry date. Inspite of above recommended storage condition, certain powdered medium may show some signs of deterioration / degradation in certain instances. This can be indicated by change in colour, change in appearance and presence of particulate matter and haziness after dissolution.

2. Preparation of concentrated medium is not recommended since free base amino acids and salt complexes having low solubility may precipitate in concentrated medium. 3. pH and sodium bicarbonate concentration of the prepared medium are critical factors affecting cell growth. This is also influenced by amount of medium and volume of culture vessel used (surface to volume ratio). For example, in large bottles, such as Roux bottles pH tends to rise perceptibly as significant volume of carbon dioxide is released. Therefore, optimal conditions of pH, sodium bicarbonate concentration, surface to volume ratio must be determined for each cell type. We recommend stringent monitoring of pH. If needed, pH can be adjusted by using sterile 1N HCl or 1N NaOH or by bubbling in carbon dioxide.

4. If required, supplements can be added to the medium prior to or after filter sterilization observing sterility precautions. Shelf life of the medium will depend on the nature of supplement added to the medium.

Disclaimer :

Revision : 04 /2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic , research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HIMEDIA