



# **Grace's Insect Medium**

With L-Glutamine and Sodium bicarbonate

# **Product Code: IML001A**

# **Product Description:**

Grace's Insect medium is a modification of the original formula developed by Wyatt. Grace (1962) improved the medium by introducing ten water soluble vitamins in Wyatt's formula. Grace successfully established four strains of cells from ovarian tissue of Australian emperor gum moth using this medium. These were the first continuous cell lines developed from insect tissue.

IML001A is Grace's Insect Medium can be used to culture cells derived from a variety of insects especially Lepidopterans and some species of Dipterans. It contains L-Glutamine and sodium bicarbonate. Originally the basal medium was supplemented with plasma from insects, the tissues of which were cultured or plasma from silkworm, *Bombyx mori*. Presently, the medium is supplemented with 5-20% fetal bovine serum. Users are advised to review the literature for recommendations regarding medium supplementation and physiological growth requirements specific for different cell lines.

## **Composition:**

Ingredients	mg/L
INORGANIC SALTS	1220.000
Calcium chloride dihydrate	1320.000
Magnesium chloride anhydrous	1068.200
Magnesium sulphate anhydrous	1356.650
Potassium chloride	2240.000
Sodium bicarbonate	350.000
Sodium phosphate monobasic	1007.000
monohydrate	
AMINO ACIDS	
DL-Serine	1100.000
Glycine	650.000
L-Alanine	225.000
L-Arginine hydrochloride	700.000
L-Asparagine monohydrate	400.000
L-Aspartic acid	350.000
L-Cystine dihydrochloride	25.000
L-Glutamic acid	600.000

L-Glutamine	600.000
L-Histidine hydrochloride monohydrate	3376.000
L-Isoleucine	50.000
L-Leucine	75.000
L-Lysine hydrochloride	625.000
L-Methionine	50.000
L-Phenylalanine	150.000
L-Proline	350.000
L-Threonine	175.000
L-Tryptophan	100.000
L-Tyrosine disodium salt dihydrate	72.000
L-Valine	100.000
ß-Alanine	200.000
VITAMINS	
Choline chloride	0.200
D-Biotin	0.010
D-Pantothenic acid (hemicalcium)	0.020
Folic acid	0.020
Niacin	0.020
Pyridoxine hydrochloride	0.020
Riboflavin	0.020
Thiamine hydrochloride	0.020
myo-Inositol	0.020
p-Amino benzoic acid (PABA)	0.020
OTHERS	
Alpha-Ketoglutaric acid	370.000
D(+) Glucose	700.000
D-Fructose	400.000
Fumaric acid free acid	55.000
L-Malic acid free acid	670.000
Succinic acid	60.000
Sucrose	26680.000

# Material required but not provided:

Fetal bovine serum (RM1112/ RM10432)

## **Quality Control:**

#### Appearance

Yellow to pale yellow colored clear solution

#### pН

6.10 - 6.30

# Osmolality in mOsm/Kg H<sub>2</sub>O 350.00 - 390.00

#### Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

#### **Cultural Response**

The growth promotion capacity of the medium is assessed qualitatively by analyzing the cells for the morphology and quantitatively by estimating the cell counts and comparing it with a control medium through minimum three subcultures.

Endotoxin Content NMT 5EU/ml

# **Storage and Shelf Life:**

Store at 2-8°C away from bright light. Shelf life is 12 months. Use before expiry date given on the product label.

#### Disclaimer :

Revision: 02/2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia<sup>™</sup> publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia<sup>™</sup> Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

