

Name of the Product : **Adult Bovine Serum
Heat Inactivated, Sterile filtered**

Code No. : **RM10913**

Section 1 : **Chemical Identification**

Code No. : RM10913

Name of the Product : Adult Bovine Serum
Heat Inactivated, Sterile filtered

Produced by : HiMedia Laboratories Pvt. Ltd.

Address : 23, Vadhani Indl. Estate, LBS Marg, Mumbai 400 086, India.

Tel. No. : 2500 0970, 2500 1607 Fax No. 022 2500 2468

Section 2 : **Adult Bovine Serum**

Section 3 : **Hazards Identification**

Hazard : Not classified as hazardous.

Section 4 : **First - Aid Measures**

No specific measures necessary.

Section 5 : **Fire Fighting Measures**

Not combustible.

Section 6 : **Accidental Release Measures**

No specific measures necessary.

Section 7 : **Handling and Storage**

Handling - Refer to Section 8

Storage - Store at -20°C

Section 8 : **Exposure Controls / Personal Protection**

Wear appropriate NIOSH/MSHA-approved respirator, chemical-resistant gloves, safety goggles, other protective clothing. Mechanical exhaust required.

Section 9 : **Physical and Chemical Properties**

Appearance : Liquid

Colour : Straw to amber

Section 10 : **Stability and Reactivity**

Stability : Materials to avoid : Strong oxidizing agents.

Hazardous decomposition products : Nature of decomposition products not known.

Hazardous polymerization : Will not occur.

Section 11 : **Toxicological Information**

Non toxic.

Section 12 : **Ecological Information**

Data not available

Section 13 : **Disposal Considerations**

No special disposal method required except that it be in accordance with current and local authority regulation.

Section 14 : **Transport Information**

UN No. : Not applicable.

Section 15 : **Regulatory Information**

Risk Phrases : Not applicable

Safety Phrases : Not applicable

Section 16 : **Other Information**

The information contained in this data sheet represents the best information currently available to us. However, no warranty is made with respect to its completeness and we assume no liability resulting from its use. The information is offered solely for user's obligation to investigate and determine the suitability of the information for their particular purpose.

