

## Filtration Funnel, Autoclavable

**Product Code: SF190A**

### Product description :

SF190A is filtration funnel made up of polysulfone which provides it chemical resistance and thermal stability upto 130°C. It is autoclavable at 121°C for 15 min at 1 bar pressure upto 100 cycles.

### Applications:

Water and product microbiology in pharmaceutical, beverages and processed food industry

### Specifications:

#### 1. Funnel

- a. Raw material : Polysulfone
- b. Capacity of funnel : 300ml
- c. Diameter : 47 mm

#### 2. Stop cork

- a. Material of construction : Silicon
- b. Size : No. 8

#### 3. Autoclavable

Autoclavable at 121°C for 15 minutes, 1 bar

### Packaging:

- Individually packed in easy-to-open bag

### Ordering information:

Code	Quantity/Pack
SF190A-INO	1No

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.