
1N Hydrochloric Acid Solution

Sterile filtered

Product Code: TCL003

Product Description :

Molecular Formula: HCl
Molecular Weight: 36.46

TCL003 is 1N hydrochloric acid solution prepared in cell culture grade water. It is used to buffer cell culture medium as per application.

Storage and Shelf Life:

Store the product at 15-30°C.
Shelf life of the product is 36 months.
Use before expiry date given on the product label.

Quality Control:

Appearance

Colorless, clear liquid.

Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

Cultural Response

No toxicity to cells.

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.