



Potassium bicarbonate Solution 7.5%

Sterile filtered

Product Code: TCL024

Product Description:

Molecular Formula: KHCO₃ Molecular Weight: 100.12

Potassium bicarbonate is used to supplement Sheilds and

Sangs M3 Insect Medium (IM004).

TCL024 is a sterile filtered solution formulated to contain 7.5% potassium bicarbonate.

Quality Control:

Appearance

Colorless, clear solution.

Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

Cultural Response

No toxicity to cells.

Endotoxin Content

NMT 5EU/ml

Storage and Shelf Life:

Store the product at 15-30°C. Shelf life of the product is 24 months. Use before expiry date given on the product label.

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMediaTM Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMediaTM Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.