



Dulbecco's Phosphate Buffered Saline 1X

With 1000mg/L D-Glucose, 36mg Sodium pyruvate per litre, Calcium and Magnesium Without Phenol red

Product Code: TL1106

Product Description:

All media used in tissue culture have a basis of a synthetic mixture of inorganic salts known as a physiological or balanced salt solution (BSS). All the physiological salt solutions have been derived from the salt solution originally described by Sydney Ringer (1885). The first balanced salt solution to be developed specifically for supporting the metabolism of mammalian cells was Tyrode's solution. Since then many modifications have been done to obtain better buffering salt solutions and to prevent calcium precipitation.

The function of a salt solution is:

- To maintain the medium within physiological pH range.
- To maintain intracellular and extra cellular osmotic balance.
- Modified with a carbohydrate, such as glucose serves as an energy source for cell metabolism.

TL1106 is Dulbecco's Phosphate Buffered Saline with 1gm glucose per litre, 36mg sodium pyruvate per litre, calcium and magnesium. It does not contain phenol red.

Composition:

Ingredients	mg/L
INORGANIC SALTS	
Calcium chloride dihydrate	133.000
Disodium hydrogen phosphate	1150.000
Magnesium chloride hexahydrate	100.000
Potassium chloride	200.000
Potassium phosphate, monobasic	200.000
Sodium chloride	8000.000
OTHERS	
D - Glucose	1000.000
Sodium pyruvate	36.000

Disclaimer :

Revision: 04/2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HIMEDIA

HiMedia Laboratories Pvt. Ltd.

Plot No. C40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane (West) 400604, Maharashtra, India.Tel No.022-69034800 Email: atc@himedialabs.com Website: www.himedialabs.com.

Quality Control:

Appearance Colorless, clear solution

pH 6.70 -7.30

Osmolality in mOsm/Kg H₂O

290.00-330.00

Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

Toxicity test Passes

Endotoxin Content NMT 1EU/ml

Storage and Shelf Life:

Store at 15-30°C away from bright light. Shelf life is 24 months. Use before expiry date given on the product label.